

### Syringe Instructions:

1. To set the syringe plunger, unlock the dose ring by rotating  $\frac{1}{4}$  turn and slide the dose ring along the plunger shaft so that the side nearest the barrel is set at the appropriate daily dose marking.
2. Rotate  $\frac{1}{4}$  turn to lock ring before dosing.
3. Make sure horse's mouth contains no feed before administration.
4. Remove syringe tip cover.
5. Insert syringe into the corner of the horse's mouth.
6. Depress the plunger until it stops at the dose ring.
7. The entire dose should be deposited on the back of the tongue or deep in the cheek pouch.



8. Horses should be observed briefly to assure no part of the dose is lost or rejected.
9. If any of the dose is lost, re-dosing is recommended.
10. Replace cap if any unused doses remain in the syringe.

